

TP53I11 Antibody




PACO23162

Product Information

Size:	Protein Background:
100ul	Induced by TP53/p53.
Reactivity:	Gene ID:
Human	TP53I11
Source:	Uniprot
Rabbit	O14683
Isotype:	Synonyms:
IgG	Tumor protein p53-inducible protein 11; p53-induced gene 11 protein; TP53I11; PIG11;
Applications:	Immunogen:
ELISA, WB, IHC	Synthesized peptide derived from internal of human TP53I11.
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, WB:1:500-1:3000, IHC:1:50-1:100	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HUVEC cells, using TP53111 antibody.

Immunohistochemistry analysis of paraffin-embedded human brain tissue using TP53111 antibody.

