

DLX3 Antibody



PACO23171

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100

Protein Background:

Likely to play a regulatory role in the development of the ventral forebrain. May play a role in craniofacial patterning and morphogenesis.

Gene ID:

DLX3

Uniprot

O60479

Synonyms:

Homeobox protein DLX-3; DLX3;

Immunogen:

Synthesized peptide derived from internal of human DLX3.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Western blot analysis of extracts from 293 cells, using DLX3 antibody.

Immunohistochemistry analysis of paraffin-embedded human heart tissue using DLX3 antibody.

