

PACO23173

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:3000,  
IF:1:100-1:500

**Protein Background:**

**Gene ID:**

SYNCRIP

**Uniprot**

O60506

**Synonyms:**

Heterogeneous nuclear ribonucleoprotein Q; hnRNP Q; hnRNP-Q; Synaptotagmin-binding; cytoplasmic RNA-interacting protein

**Immunogen:**


Synthesized peptide derived from Internal of human hnRNP Q.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from HUVEC cells, Jurkat cells and 293 cells, using hnRNP Q antibody.

Immunofluorescence analysis of HepG2 cells, using hnRNP Q antibody.

