

PTTG1IP Antibody



PACO23232

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100

Protein Background:

May facilitate PTTG1 nuclear translocation.

Gene ID:

PTTG1IP

Uniprot

P53801

Synonyms:

C21orf1; PBF; PTTG; PTTG-binding factor; pituitary tumor-transforming gene 1 protein-interacting protein


Immunogen:

Synthesized peptide derived from C-terminal of human PTTG.


Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from RAW264.7 cells, using PTTG antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using PTTG antibody.