

SOX8/SOX9/SOX17/SOX18 Antibody



PACO23234

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100

Protein Background:

May play a role in central nervous system, limb and facial development. May be involved in male sex determination. Binds the consensus motif 5'-[AT][AT]CAA[AT]G-3' By similarity.

Gene ID:

SOX8/SOX9/SOX17/SOX18

Uniprot

P57073/P48436/Q9H6I2/P35713

Synonyms:

Transcription factor SOX-8/9/17/18; SOX8/9/17/18;


Immunogen:

Synthesized peptide derived from C-terminal of human SOX8/9/17/18.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from COS7 cells, untreated or treated with serum (10%, 30mins), using SOX8/9/17/18 antibody.

Immunohistochemistry analysis of paraffin-embedded human brain tissue using SOX8/9/17/18 antibody.

