

KSR2 Antibody



PACO23257

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IF:1:100-1:500

Protein Background:**Gene ID:**

KSR2

Uniprot

Q6VAB6

Synonyms:

hKSR2; Kinase suppressor of Ras 2; KSR2;


Immunogen:

Synthesized peptide derived from internal of human KSR2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from COS7, 293 and 3T3 cells, using KSR2 antibody.

Immunofluorescence analysis of HeLa cells, using KSR2 antibody.

