

# TP53AIP1 Antibody



PACO23348

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100,  
IF:1:100-1:500

**Protein Background:**

May play an important role in mediating p53/TP53-dependent apoptosis.

**Gene ID:**

TP53AIP1

**Uniprot**

Q9HCN2

**Synonyms:**

p53-regulated apoptosis-inducing protein 1; P53aip1; TP53AIP1;

**Immunogen:**

Synthesized peptide derived from internal of human TP53AIP1.


**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunohistochemistry analysis of paraffin-embedded human thymus gland tissue using TPIP1 antibody.



Immunofluorescence analysis of A549 cells, using TPIP1 antibody.

