

EPHB1/EPHB2/EPHB3 Antibody



PACO23432

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100,
IF:1:100-1:500

Protein Background:**Gene ID:**

EPHB1/EPHB2/EPHB3

Uniprot

P54762/P29323/P54753

Synonyms:

EC 2.7.10.1; ELK; EPB1; EPHT2; EPTH2

Immunogen:


Synthesized peptide derived from internal of human EPHB1/2/3.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue using EPHB1/2/3 antibody.



Immunofluorescence analysis of NIH/3T3 cells, using EPHB1/2/3 antibody.

