

HOXA11/HOXD11 Antibody



PACO23443

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IF:1:100-1:500

Protein Background:

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Gene ID:

HOXA11/HOXD11

Uniprot

P31270/P31277

Synonyms:

Homeobox protein Hox-A11; Hox-1l; HOXA11; HOX11;


Immunogen:

Synthesized peptide derived from internal of human HOXA11/D11.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from RAW264.7 cells, using HOXA11/D11 antibody.

Immunofluorescence analysis of HUVEC cells, using HOXA11/D11 antibody.

