

USP50 Antibody



PACO23444

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100

Protein Background:

Has no peptidase activity.

Gene ID:

USP50

Uniprot

Q70EL3

Synonyms:

Inactive ubiquitin carboxyl-terminal hydrolase 50; Inactive ubiquitin-specific peptidase 50;


Immunogen:

Synthesized peptide derived from C-terminal of human USP50.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from A549 cells, using USP50 antibody.

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using USP50 antibody.

