

USP13 Antibody



PACO23492

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100

Protein Background:**Gene ID:**

USP13

Uniprot

Q92995

Synonyms:

Ubiquitin carboxyl-terminal hydrolase 13; EC 3.1.2.15; Ubiquitin thioesterase 13;
Ubiquitin-specific-processing protease 13; Deubiquitinating enzyme 13


Immunogen:

Synthesized peptide derived from internal of human USP13.


Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM
NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from mouse brain cells, using USP13 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using USP13 antibody.