

THOC5 Antibody



PACO23494

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100

Protein Background:**Gene ID:**

THOC5

Uniprot

Q13769

Synonyms:

THO complex subunit 5 homolog; NF2/meningioma region protein pK1.3; Placental protein 39.2; THOC5; C22orf19


Immunogen:

Synthesized peptide derived from internal of human THOC5.


Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from JK cells and HepG2 cells, using THOC5 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using THOC5 antibody.