

ATF2 (Ab-71 or 53) Antibody



PACO23540

Product Information

Size:

100ul

Protein Background:**Reactivity:**

Human, Mouse, Rat

Gene ID:

ATF2

Source:

Rabbit

Uniprot

P15336

Isotype:

IgG

Synonyms:

CREB2; CREBP1;

Applications:

ELISA, WB, IHC, IF

Immunogen:

Peptide sequence around aa.69~73 or 51~55 (T-P-T-P-T) derived from Human ATF2.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000,
IHC:1:50-1:200, IF:1:100-1:200


Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4,
150mM NaCl, 0.02% sodium azide and 50% glycerol.


Product Images



Western blot analysis of extracts from HT29 and K562 cells using ATF2(Ab-71 or 53) Antibody.



Immunofluorescence staining of methanol-fixed HeLa cells using ATF2(Ab-71 or 53) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using ATF2(Ab-71 or 53) Antibody(left) or the same antibody preincubated with blocking peptide(right).