

TP53 (Ab-15) Antibody



PACO23552

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000,
IHC:1:50-1:200, IF:1:100-1:200

Protein Background:**Gene ID:**

TP53

Uniprot

P04637

Synonyms:

Tumor suppressor p53; Phosphoprotein p53; Antigen NY-CO-13; TP53;


Immunogen:

Peptide sequence around aa.13~17 (P-L-S-Q-E) derived from Human p53.


Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4,
150mM NaCl, 0.02% sodium azide and 50% glycerol.


Product Images



Western blot analysis of extracts from MDA and JK cells using p53(Ab-15) Antibody.



Immunofluorescence staining of methanol-fixed HeLa cells using p53(Ab-15) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using p53(Ab-15) Antibody(left) or the same antibody preincubated with blocking peptide(right).