

EIF4G1 (Ab-1232) Antibody



PACO23594

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000,
IHC:1:50-1:200, IF:1:100-1:200

Protein Background:**Gene ID:**

EIF4G1

Uniprot

Q04637


Synonyms:**Immunogen:**

Peptide sequence around aa.1230~1234 (P-V-S-P-L) derived from Human eIF4G.


Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Product Images



Western blot analysis of extracts from HepG2 and HT29 cells using eIF4G(Ab-1232) Antibody.



Immunofluorescence staining of methanol-fixed HeLa cells using eIF4G(Ab-1232) Antibody.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using eIF4G(Ab-1232) Antibody(left) or the same antibody preincubated with blocking peptide(right).