

FGFR3 Antibody



PACO23656

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100

Protein Background:**Gene ID:**

FGFR3

Uniprot

P22607

Synonyms:

2.7.1.112; CEK2; HSFGR3EX; fibroblast growth factor receptor 3 isoform 1 precursor;
fibroblast growth factor receptor 3 isoform 2 precursor

Immunogen:


Synthesized peptide derived from human FGFR3.

Storage:


Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM
NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using FGFR3 antibody.



Western blot analysis of extracts from 293 cells (Lane 2) and HepG2 cells (Lane 3), using FGFR3 antibody. The lane on the left is treated with synthesized peptide.



Western blot analysis of extracts from LOVO cells, using FGFR3 antibody.

