

TRA2A Antibody



PACO23749

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100, IF:1:100-1:500

Protein Background:

Sequence-specific RNA-binding protein which participates in the control of pre-mRNA splicing.

Gene ID:

TRA2A

Uniprot

Q13595

Synonyms:

Transformer-2 protein homolog; TRA-2 alpha;


Immunogen:

Synthesized peptide derived from internal of human TRA-2 alpha.


Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Product Images



Western blot analysis of extracts from HUVEC cells, using TRA-2α antibody.



Immunofluorescence analysis of HeLa cells, using TRA-2 alpha antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using TRA-2 alpha antibody.