

HIPK4 Antibody



PACO23758

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100, IF:1:100-1:500

Protein Background:

Protein kinase that phosphorylates human TP53 at Ser-9, and thus induces TP53 repression of BIRC5 promoter By similarity. May act as a corepressor of transcription factors Potential.

Gene ID:

HIPK4

Uniprot

Q8NE63

Synonyms:

Homeodomain-interacting protein kinase 4; EC 2.7.11.1; HIPK4;


Immunogen:

Synthesized peptide derived from internal of human HIPK4.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HUVEC cells, using HIPK4 antibody.

Immunofluorescence analysis of A549 cells, using HIPK4 antibody.



Immunohistochemistry analysis of paraffin-embedded human testis tissue using HIPK4 antibody.

