

MPC1 Antibody



PACO23798

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100, IF:1:100-1:500

Protein Background:

Mediates the uptake of pyruvate into mitochondria.

Gene ID:

MPC1

Uniprot

Q9Y5U8

Synonyms:

BR44L; brain protein 44-like;


Immunogen:

Synthesized peptide derived from internal of human BRP44L.


Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Product Images



Western blot analysis of extracts from COLO cells and HeLa cells, using BRP44L antibody.



Immunofluorescence analysis of A549 cells, using BRP44L antibody.



Immunohistochemistry analysis of paraffin-embedded human heart tissue using BRP44L antibody.