

Phospho-CTNNB1 (Thr41/Ser45) Antibody



PACO23855

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

CTNNB1

Uniprot

P35222

Synonyms:

CTNNB1; CATNB; CTNB1; CTNNB;


Immunogen:

Peptide sequence around phosphorylation site of threonine 41/serine 45 (A-T-T(p)-T-A-P-S(p)-L-S) derived from Human b-Catenin.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from SW626 cells using β -Catenin(Phospho-Thr41/Ser45) Antibody(Lane 2) and the same antibody preincubated with blocking peptide(Lane1).