

Phospho-HDAC4 (Ser632) Antibody



PACO23864

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

HDAC4

Uniprot

P56524

Synonyms:

HD4


Immunogen:

Peptide sequence around phosphorylation site of serine 632 (A-Q-S(p)-S-P) derived from Human HDAC4.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 cells untreated(lane 1) or treated with EGF(lane 2) using HDAC4(Phospho-Ser632) Antibody.