

Phospho-EIF2AK2 (Thr451) Antibody



PACO23883

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:**Gene ID:**

EIF2AK2

Uniprot

P19525

Synonyms:

ADRB2; E2AK2; EIF2AK2; EIF2aK; PRKR


Immunogen:

Peptide sequence around phosphorylation site of threonine 451 (K-G-T(p)-L-R) derived from Human PKR.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded human colon carcinoma tissue using PKR(Phospho-Thr451) Antibody(left) or the same antibody preincubated with blocking peptide(right).