

# Phospho-NTRK2 (Tyr705) Antibody



PACO23891

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:1000

**Protein Background:****Gene ID:**

NTRK2

**Uniprot**

Q16620

**Synonyms:**

BDNF/NT-3 growth factors receptor precursor; EC 2.7.10.1; GP145-TrkB; GP145-TrkB/GP95-TrkB; NTRK2

**Immunogen:**


Peptide sequence around phosphorylation site of tyrosine 705 (T-D-Y P-Y-R) derived from Human TrkB.

**Storage:**

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from mouse brain tissue using TrkB(Phospho-Tyr705) Antibody and the same antibody preincubated with blocking peptide.