

Phospho-PIK3R1/PIK3R3 (Tyr467/199) Antibody



PACO23894

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

PIK3R1/PIK3R3

Uniprot

P27986/Q92569

Synonyms:

p85, AGM7, p85-ALPHA, p55, p55-GAMMA


Immunogen:

Peptide sequence around phosphorylation site of tyrosine 467 (L-Y(p)-E-E-Y) derived from Human PI3 Kinase p85/p55.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 cells untreated or treated with sorbitol using PI3 Kinase p85/p55 (phospho-Tyr467/199) Antibody.