

Phospho-ACACA (Ser79) Antibody



PACO23917

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ACACA

Uniprot

Q13085


Synonyms:

ACC, ACAC, ACC1, ACCA,

Immunogen:**Storage:**

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 cells untreated or treated with AICAR using Acetyl-CoA Carboxylase (Phospho-Ser79) Antibody.