

Phospho-TP53 (Ser20) Antibody



PACO23940

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

TP53

Uniprot

P04637

Synonyms:

P53, BCC7, LFS1, TRP53,


Immunogen:

Peptide sequence around phosphorylation site of Serine 20(T-F-S(p)-D-L) derived from Human p53.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Mouse brain tissue using p53 (Phospho-Ser20) Antibody. The lane on the right is treated with the antigen-specific peptide.