

Phospho-RPS6KB1 (Ser371) Antibody



PACO23942

Product Information

Size:

100ul

Protein Background:**Reactivity:**

Human, Mouse, Rat

Gene ID:

RPS6KB1

Source:

Rabbit

Uniprot

P23443

Isotype:

IgG

Synonyms:

S6K, PS6K, S6K1, STK14A, p70-S6K

Applications:

ELISA, WB

Immunogen:

Peptide sequence around serine 371 (V-D-S(p)-P-D) derived from Human p70 S6 Kinase.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Mouse kidney tissue using p70 S6 Kinase (Ser371) Antibody. The lane on the left is treated with the antigen-specific peptide.