

**Phospho-MAPK3/MAPK1 (Thr202/Tyr204)
Antibody
PACO23995**



Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:

Gene ID:

MAPK3/MAPK1

Uniprot

P27361/P28482

Synonyms:

ERK1, ERT2, ERK-1, PRKM3, P44ERK1


Immunogen:

Peptide sequence around phosphorylation site of Threonine 202/Tyrosine 204(F-L-T(p)-E-Y(p)-V-A)Sderived from Human ERK1/2.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 3T3 cells treated with IFN α using ERK1/2 (Phospho-Thr202/Tyr204) Antibody. The lane on the left is treated with the antigen-specific peptide.