

Phospho-ABL1 (Tyr204) Antibody



PACO23998

Product Information

Size:

100ul

Protein Background:**Reactivity:**

Human, Mouse, Rat

Gene ID:

ABL1

Source:

Rabbit

Uniprot

P00519

Isotype:

IgG

Synonyms:

p150; c-ABL; EC 2.7.10.2;

Applications:

ELISA, WB

Immunogen:

Peptide sequence around phosphorylation site of tyrosine 204 (K-L-Y(p)-V-S) derived from Human c-Abl.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from COS7 cells treated with Adriamycin using Abl (Phospho-Tyr204) antibody. The lane on the right is treated with the antigen-specific peptide.