

# Phospho-CD19 (Tyr531) Antibody



PACO24006

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:1000

**Protein Background:****Gene ID:**

CD19

**Uniprot**

P15391

**Synonyms:**

B-lymphocyte antigen CD19; B-lymphocyte surface antigen B4; CD19; CD19 molecule; MGC12802

**Immunogen:**


Peptide sequence around phosphorylation site of tyrosine 531 (D-S-Y(p)-E-N) derived from Human CD19.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from COS7 cells treated with Serum using Phospho-CD19 (Tyr531) antibody. The lane on the right is treated with the antigen-specific peptide.