

Phospho-PDGFRα (Tyr849) Antibody



PACO24008

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

PDGFRA

Uniprot

P16234

Synonyms:

CD140a; PDGFR-alpha; PGFRA; kinase PDGFR-alpha;


Immunogen:

Peptide sequence around phosphorylation site of Tyrosine849(S-N-Y(p)-V-S) derived from Human PDGFRA.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 tissue using PDGFRα (Phospho-Tyr849) antibody. The lane on the right is treated with the antigen-specific peptide.