

Phospho-FGFR1 (Tyr766) Antibody



PACO24032

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

FGFR1

Uniprot

P11362

Synonyms:

FGFBR; FGFR-1; FLG; FLT2; c-fgr


Immunogen:

Peptide sequence around phosphorylation site of tyrosine 766 (Q-E-Y(p)-L-D) derived from FGFR1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 cells treated with EGF using Phospho-FGFR1 (Tyr766) antibody. The lane on the right is treated with the antigen-specific peptide.