

# Phospho-NTRK1 (Tyr791) Antibody



PACO24083

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Protein Background:****Gene ID:**

NTRK1

**Uniprot**

P04629

**Synonyms:**

TRKA; TRK; NTRK1; p140-TrkA;

**Immunogen:**


Peptide sequence around phosphorylation site of tyrosine 791 (P-V-Y(p)-L-D) derived from Human Trk A.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Immunohistochemical analysis of paraffin-embedded human brain tissue, using Trk A (phospho-Tyr791) antibody (left) or the same antibody preincubated with blocking peptide (right).