

# Phospho-IL2RB (Tyr364) Antibody



PACO24096

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:1000

**Protein Background:****Gene ID:**

IL2RB

**Uniprot**

P14784

**Synonyms:**

CD122 antigen; P70-75; IL-2 receptor; IL-2R beta;

**Immunogen:**


Peptide sequence around phosphorylation site of tyrosine 364(Q-G-Y(p)-F-F) derived from Human IL-2R &beta; /CD122.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from COS7 cells using IL-2R beta /CD122 (Phospho-Tyr364) Antibody. The lane on the right is treated with the antigen-specific peptide.