

Phospho-MAP2K7 (Thr275) Antibody



PACO24116

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

MAP2K7

Uniprot

O14733

Synonyms:

PRKMK7; JNKK2; MKK7; MAP2K7; MP2K7


Immunogen:

Peptide sequence around phosphorylation site of threonine 275(A-K-T(p)-R-S) derived from Human MAP2K7.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from cos-7 cells (Lane 2) and 3T3 cells (Lane 3), using MAP2K7 (Phospho-Thr275) Antibody. The lane on the left is treated with antigen-specific peptide.