

Phospho-INSR (Thr1375) Antibody



PACO24188

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

INSR

Uniprot

P06213

Synonyms:

IR; CD220 antigen; kinase InsR;


Immunogen:

Peptide sequence around phosphorylation site of threonine 1375(I-L-T(p)-L-P) derived from Human IR.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 cells using INSR (Phospho-Thr1375) Antibody. The lane on the right is treated with the antigen-specific peptide.