

Phospho-PTK2 (Tyr576/Tyr577) Antibody



PACO24345

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000,
IF:1:100-1:200

Protein Background:**Gene ID:**

PTK2

Uniprot

Q05397

Synonyms:

FADK 1; FAK1; PTK2


Immunogen:

Peptide sequence around phosphorylation site of tyrosine 576/tyrosine 577 (S-T-Y(p)-Y(p)-K-A) derived from Human FAK.


Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HT29 cells untreated or treated with PMA using FAK(phospho-Tyr576/Tyr577) Antibody.



Immunofluorescence staining of methanol-fixed Hela cells using FAK(phospho-Tyr576/Tyr577) Antibody.