

Phospho-ABL1 (Tyr412) Antibody



PACO24355

Product Information

Size:

100ul

Protein Background:**Reactivity:**

Human, Mouse, Rat

Gene ID:

ABL1

Source:

Rabbit

Uniprot

P00519

Isotype:

IgG

Synonyms:

ABL, JTK7, p150, c-ABL,

Applications:

ELISA, IHC, IF


Immunogen:**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100,
IF:1:100-1:200


Storage:

Purified Rabbit polyclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Product Images



Immunofluorescence staining of methanol-fixed HeLa cells using c-Abl (phospho-Tyr412) antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using c-Abl (phospho-Tyr412) antibody.