

Phospho-CARM1 (Ser228) Antibody



PACO24363

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000,
IF:1:100-1:200

Protein Background:**Gene ID:**

CARM1

Uniprot

A6NN38

Synonyms:

PRMT4


Immunogen:

Peptide sequence around phosphorylation site of serine 228(V-K-S(p)-N-N) derived from Human CARM1.


Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from A431 cells untreated or treated with EGF (200ng/ml, 5min), using CARM1 (Ab-228) antibody (Line 1 and 2) and CARM1 (Phospho-Ser228) antibody (Line 3 and 4).



Immunofluorescence staining of methanol-fixed HeLa cells using CARM1 (Phospho-Ser228) antibody.