

Phospho-TP53 (Ser33) Antibody



PACO24456

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000,
IF:1:100-1:200

Protein Background:**Gene ID:**

TP53

Uniprot

P04637

Synonyms:

P53, BCC7, LFS1, TRP53


Immunogen:

Peptide sequence around phosphorylation site of serine 33 (V-L-S(p)-P-L) derived from Human p53.


Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Product Images



Western blot analysis of extracts from HT29 cells using p53(Phospho-Ser33) Antibody and the same antibody preincubated with blocking peptide.



Western blot analysis of extracts from HT29 cells, treated with serum or calf intestinal phosphatase (CIP), using p53 (Phospho-Ser33) Antibody.



Immunofluorescence staining of methanol-fixed HeLa cells using p53(Phospho-Ser33) Antibody.