

PACO24610

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

WB

Recommended dilutions:

WB:1:500-1:3000

Protein Background:

This is a copper-containing oxidase that functions in the formation of pigments such as melanins and other polyphenolic compounds. Catalyzes the rate-limiting conversions of tyrosine to DOPA DOPA to DOPA-quinone and possibly 56-dihydroxyindole to indole-56 quinone.

Gene ID:

TYR

Uniprot

P14679

Synonyms:

TYR; EC 1.14.18.1; Monophenol monooxygenase; Tumor rejection antigen AB; SK29-AB

Immunogen:

The antibody detects endogenous levels of total tyrosinase protein.

Storage:

pH 7.4 150mM NaCl 0.02% sodium azide and 50% glycerol. al

Product Images

N/A

N/A