

PACO24628

---

## Product Information

**Size:**

50ug

**Reactivity:**

Chicken

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Serum albumin, the main protein of plasma, has a good binding capacity for water, Ca(2+), Na(+), K(+), fatty acid, , hormones, bilirubin and drugs. Its main function is the regulation of the colloidal osmotic pressure of blood. Major zinc transporter in plasma, typically binds about 80% of all plasma zinc.

**Gene ID:**

ALB

**Uniprot**

P19121

**Synonyms:**

Serum albumin (Alpha-livetin) (allergen Gal d 5), ALB

**Immunogen:**

Serum Albumin purified from Chicken serum (Native Protein).

**Storage:**

Preservative: 0.1mol/L NaCl, 15mmol/L NaN<sub>3</sub>, PH7.2

## Product Images

---

N/A

N/A