

RPL11 Antibody, HRP conjugated



PACO24653

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RPL11

Uniprot

P62913

Synonyms:

60S ribosomal protein L11 (CLL-associated antigen KW-12) (Large ribosomal subunit protein uL5), RPL11

Immunogen:

Recombinant Human 60S ribosomal protein L11 protein (3-178AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A