

PACO24656

Product Information

Size:

50ug

Reactivity:

Rabbit

Source:

Guinea pig

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Serum albumin, the main protein of plasma, has a good binding capacity for water, Ca(2+), Na(+), K(+), fatty acid, , hormones, bilirubin and drugs. Its main function is the regulation of the colloidal osmotic pressure of blood. Major zinc transporter in plasma, typically binds about 80% of all plasma zinc.

Gene ID:

ALB

Uniprot

P49065

Synonyms:

Serum albumin, ALB

Immunogen:

Serum Albumin purified from Rabbit serum (Native Protein).

Storage:

Preservative: 0.1mol/L NaCl, 15mmol/L NaN₃, PH7.2

Product Images

N/A

N/A