

# A4GNT Antibody, HRP conjugated



PACO24669

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Necessary for the synthesis of type III mucin. Catalyzes the transfer of N-acetylglucosamine (GlcNAc) to core 2 branched O-glycans.

**Gene ID:**

A4GNT

**Uniprot**

Q9UNA3

**Synonyms:**

Alpha-1,4-N-acetylglucosaminyltransferase (Alpha4GnT) (EC 2.4.1. -), A4GNT

**Immunogen:**

Recombinant Human Alpha-1,4-N-acetylglucosaminyltransferase protein (26-340AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A