

PACO24684

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Mediates ATP-dependent export of organic anions and drugs from the cytoplasm. Hydrolyzes ATP with low efficiency. Not capable of conferring drug resistance. Mediates the translocation of phosphatidylcholine across the canalicular membrane of the hepatocyte.

Gene ID:

Abcb4

Uniprot

P21440

Synonyms:

Phosphatidylcholine translocator ABCB4 (ATP-binding cassette sub-family B member 4) (Multidrug resistance protein 2) (Multidrug resistance protein 3) (EC 3.6.3.44) (P-glycoprotein 2) (P-glycoprotein 3), Abcb4, Mdr2 Pgy-2 Pgy2

Immunogen:

Recombinant Mouse Phosphatidylcholine translocator ABCB4 protein (352-708AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A