

ACAA1 Antibody



PACO24700

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

ACAA1

Uniprot

P09110

Synonyms:

3-ketoacyl-CoA thiolase, peroxisomal (EC 2.3.1.16) (Acetyl-CoA acyltransferase) (Beta-ketothiolase) (Peroxisomal 3-oxoacyl-CoA thiolase), ACAA1, ACAA PTHIO


Immunogen:

Recombinant Human 3-ketoacyl-CoA thiolase, peroxisomal protein (27-331AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO24700 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO24700 at dilution of 1:100.