

PACO24716

Product Information

Size:

50ug

Reactivity:

Bovine

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

erminates signal transduction at the neuromuscular junction by rapid hydrolysis of the acetylcholine released into the synaptic cleft.

Gene ID:

ACHE

Uniprot

P23795

Synonyms:

Acetylcholinesterase (AChE) (EC 3.1.1.7), ACHE

Immunogen:

Recombinant Bovine Acetylcholinesterase protein (31-613AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A