

ACTRT1 Antibody



PACO24739

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

mitophagy in response to mitochondrial depolarization, positive regulation of defense response to virus by host, xenophagy.

Gene ID:

ACTRT1

Uniprot

Q8TDG2

Synonyms:

Actin-related protein T1 (ARP-T1), ACTRT1, ARPT1


Immunogen:

Recombinant Human Actin-related protein T1 protein (201-376AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human colon tissue using PACO24739 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO24739 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO24739 at dilution of 1:100.