

ACTRT2 Antibody



PACO24743

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

ACTRT2

Uniprot

Q8TDY3

Synonyms:

Actin-related protein T2 (ARP-T2) (Actin-related protein M2), ACTRT2, ARPM2


Immunogen:

Recombinant Human Actin-related protein T2 protein (1-376AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: ACTRT2 antibody at 6 μ g/ml

Lane 1: HepG2 whole cell lysate

Lane 2: A431 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 42 kDa

Observed band size: 42 kDa

.